



Empore™

High Performance Extraction Disk Plates

PRODUCT ———
—— Information

Instructions for Use


Product Description

Empore™ Extraction Disk Plates are designed for high throughput solid phase extraction (SPE). 96 samples can be processed within a standard 8 row by 12 column microtiter plate format. One disk plate can replace four separate runs on a conventional SPE manifold handling 24 individual cartridges per run. The 96 well format is ideal for sample preparation prior to LC/MS/MS or other high throughput analytical techniques.

The plate is molded from a polypropylene resin suggested for medical and pharmaceutical applications. An Empore extraction disk is secured in place at the bottom of each well with a sealing ring. A patented prefilter is placed above the Empore disk. This prefilter aids in preventing particulates and macromolecules from reaching the underlying membrane and improves the flow of biological samples, such as serum and plasma, through the plate.

The prefilter is composed of polypropylene microfiber layers of graded densities. Three different densities are used, with the coarsest one on top and the finest at the bottom. The top two microfiber layers are individual layers of material. The third microfiber layer, having the smallest effective pore size, is on the bottom of the prefilter and contains five individual layers of material. A porous polypropylene support membrane comprises the final layer.

Empore™ Extraction Disk Plate



Product Characteristics

Each well of the Empore™ Extraction Disk Plate has an effective membrane diameter of 5.5 mm. Effective diameter is the diameter of membrane available for use during sample processing. This diameter is smaller than the actual well diameter due to the dimensions of the rings that seal the membrane and prefilter into place. The reservoir volume of the standard plate is 1.2 mL and the volume of the deep well plate is 2.5 mL. If sample or reagent volumes exceed the volume of the well, multiple aliquots of solution may be used.

Standard density (SD) membranes are composed of chromatographic particles commonly referred to as from 40-60 μm in size (actual mean size is about 55 μm). The standard density membrane has been optimized for improved flow rates for samples processed in most bioanalytical applications.

Typical Empore™ Extraction Disk Plate Specifications

Typical Empore™ Extraction Disk Plate Specifications	
Membrane diameter	5.5 mm
Well volume	1.2 mL and 2.5 mL
Membrane thickness	0.75 mm
Membrane type	Standard density (SD)
Prefilter composition	Graded density polypropylene
Bed volume	18 μL
Bonded silica sorbent mass (C2, C8, C18 and MPC Mixed Phase Cation)	10 mg (nominal)
Universal Resin sorbent mass	5 mg (nominal)
Mean particle size	55 μm (bonded silicas), 44 μm (Universal Resin)
Membrane composition	<ul style="list-style-type: none">• Bonded silicas: nominally 93% sorbent, 7% PTFE (w/w)• Universal Resin: nominally 90% sorbent, 10% PTFE (w/w)

Empore™ Disk Technology

Empore™ solid phase extraction disks are produced by trapping sorbent particles within an inert matrix of polytetrafluoroethylene (PTFE). The resulting particle-loaded membrane yields a denser, more uniform extraction bed than can be achieved with traditional loosely packed SPE particles. The result is improved mass transfer kinetics with consistent performance in solid phase extraction methods.

The dense particle packing and uniform distribution within Empore disks offer outstanding sample preparation efficiency and reproducibility of results. Since the diffusion distance between particles is minimized, adsorption is more efficient, and extraction can be accomplished using low sorbent mass. The following performance gains can be realized:

- Reduced solvent volumes
- Small elution volumes
- Reduced time for eluate evaporation
- Potential elimination of eluate evaporation
- Higher throughput
- Channeling effects eliminated
- Excellent reproducibility/low CVs



Cross section of a standard density (SD) Empore™ disk (40-60 μm particle size)

Plate Processing Equipment

The dense particle packing of the disk prevents samples from flowing by gravity through the membrane. Vacuum, positive pressure, and centrifugation can all be used to process liquids through the membrane.

Vacuum Manifold

The Empore extraction disk plate differs in design from plates offered by other manufacturers and requires a specific vacuum manifold that is supplied by 3M (Cat. 610 Vacuum Manifold). Manifolds supplied by different plate manufacturers are specific to the design of their respective plates and cannot generally be used with plates from other manufacturers. In some cases, the use of spacers to raise/lower the height of a collection device inside the manifold may overcome some incompatibility, but must specifically be determined for each combination of manufacturer's manifold and the collection device.

The Empore™ Vacuum Manifold (Cat. 610) is illustrated at the right. A waste tray is shown that is used for liquid collection during the conditioning, loading and washing steps of the extraction process.

During the elution step, several types of collection devices may be used. These include deep well plates, shallow well plates, or racks of collection tubes. Please refer to the "Collection Device Options" section for more information about collection plates.

Plates can be processed manually on a vacuum manifold using multichannel pipettors (8- or 12-channel) for liquid transfers. The basic items needed include:

- A rack of pipette tips.
- Several disposable reservoir trays to hold each liquid (e.g., conditioning, wash and elution solvents).
- A collection device (plate or microtube format) for the eluate.
- A placeholder, such as an 8-well cover strip, may also be useful.

If samples are not contained in the 8 x 12 microtiter format before use, individual transfers will be necessary from vials or tubes into respective positions within the disk plate.



Plate Processing Equipment (cont.)

Centrifugation

Centrifugation is a very effective method for processing liquids through the plate. The high forces attainable with centrifugation are advantageous to displace liquids after each step (e.g., to fully remove aqueous wash before the elution step). The Empore™ extraction disk plate is nested within a collection device, such as a deep well plate, or rack of microtubes, and this combination is inserted into a microplate carrier that is connected to the rotor. Centrifuges are available that process from two to four microplates per run.

Centrifugation times and speeds should be optimized for individual samples and analytes. Sample loading and eluting should be evaluated at different G forces to check the impact on analyte recovery. Some general guidelines for getting started are found below:

Sample	Sample Extraction	Elution
Buffer	100 G/2 minutes	100 G/2 minutes
Serum diluted 1:3	100 G/2 minutes	100 G/2 minutes
Serum diluted 1:1	250 G/2 minutes	100 G/2 minutes
Undiluted serum	250 G/2 minutes	100 G/2 minutes



Centrifugation is an effective method to process liquids through the Empore™ extraction disk plate

Liquid Handling Robots

The throughput gains provided by the Empore extraction disk plate attain an even higher level of performance when the extraction is automated, thereby eliminating the possibility for human pipetting error and/or procedural error. Commercially available liquid handling workstations provide a variety of automation solutions since many are designed to work with the microtiter plate format.

Vacuum Recommendations

Empore extraction disk plates generally require from 10-20 inches Hg (0.34 to 0.68 bar) vacuum to process biological fluids such as plasma and serum through all 96 wells of the plate. If the sample matrix is relatively clean and of a small volume, lower vacuum may be used. However, a general guide is to open the vacuum source to about 17 inches Hg (0.58 bar) for all steps, with the following exceptions:

1. During method optimization, try loading the sample matrix at both a low vacuum (5-7 inches Hg; 0.17 to 0.24 bar) and a high vacuum (about 17-20 inches Hg; 0.58 to 0.68 bar). Compare recoveries at both levels of vacuum and select the vacuum setting that optimizes recoveries. If an analyte has a low affinity for the sorbent, it may need to pass through the sorbent bed more slowly during the loading step for sufficient sorption to occur.
2. Lower vacuum during the elution step will yield a slower flow rate, which is important to prevent splashing in shallow well collection plates.
3. Lower vacuum/slower flow rate may also be beneficial when eluting from MPC sorbent, as more time may be needed to disrupt ionic interactions.



Maintaining Uniform Flow


To maintain uniform flow between wells and avoid loss of vacuum caused by some wells emptying faster than others, several approaches should be considered.

1. During sample loading, inconsistent flow is generally caused by differences in sample viscosities between wells within a 96 well plate. If possible, avoid loading buffer and plasma samples within the same plate. Viscosity differences can also be addressed by diluting serum and plasma samples. Depending on sample size and well volume within the plate, plasma or serum samples may be diluted up to 1:4 with buffer.
2. During the wash and elution steps, it is important to avoid precipitating proteins that can reduce flow through the disk. The most effective method of eliminating proteins is to follow the sample load step immediately with a 100% water wash before washing with solutions containing a small percentage of organic. Please refer to the Empore™ High Performance Extraction Disk Plates Method Optimization Guide for more detailed information.
3. If the above techniques fail to resolve the problem, centrifugation of plates and the use of positive pressure may be considered as alternative processing approaches.

Sealing Unused Wells

Empore™ Cat. 660 Sealing Tape Pads are suggested to cover unused wells during an extraction assay. Use a clean sheet each time for maximum sealing ability. Do not try to reuse a sheet. *Note that these sheets are not recommended for sealing collection plates before injection.* After a partial plate is used, an additional sheet can be employed to cover the used portion so that it is not mistakenly used again.


Extraction disk plates are designed as a single use consumable item and 3M does not warrant their multiple use.



Sealing Unused Wells
with Cat. 660 Sealing Tape Pad

Collection Device Options


The illustration below represents some typical collection device options for use with the 96-well extraction disk plate. Cross contamination can be prevented by fitting the collars of the Empore™ Extraction Disk Plate inside the well of the collection plate. Polypropylene shallow well plates are used for very small elution volumes and are available in either a round U-bottom or a tapered V-bottom configuration. Microtubes are individual polypropylene tubes arranged in a rack conforming to the 8 x 12 microtiter array. Individual assay conditions, elution volumes and subsequent eluate handling procedures define the selection of a collection device.




Elution Using Glass Inserts

Other collection device options, such as smaller volume microtubes and smaller volume deep well plates can be used. Microtubes are another option and are available in strips of eight connected tubes as well as individual tubes.

If glass is preferred as a collection device instead of polypropylene, deep and shallow well plates can be used with glass inserts, either capped (PTFE septum) or uncapped.






The Empore™ Cat. 610 Vacuum Manifold accepts deep and shallow well collection plates and microtubes. Spacers supplied with the manifold may be used to raise the collection device to the optimal distance below the tips of the disk plate. The illustrations below show that two spacers are used with a deep well plate, no spacers are necessary for the microtubes, and a shallow well plate uses an empty microtube rack as a spacer.



Optimal Heights of Plates with Various Collection Devices

The optimal distances between disk plate tip and collection plate are shown in the accompanying illustrations. Using two spacers, the collar of the disk plate fits into the wells of the deep well plate before the vacuum is applied. After vacuum is turned on, the plate will lower further. These conditions help to avoid cross contamination and ensure eluate integrity.

A shallow well plate does not hold much volume, so it is important to use this plate with caution to avoid overfilling the wells. The tips of the disk plate should fit just into the top of the well. Once vacuum is turned on, the tips will lower themselves a bit more. *The important point to remember is that the tip of the plate must not touch any liquid collected in the well!* Bubble formation can cause liquid splashing out of the well and into an adjacent well. Shallow well plates work best with eluate volumes of 150 μL or less.



Autoinjectors are commercially available that accept deep and/or shallow well plates. After covering the plate with an appropriate material for sealing, it can be loaded into the injector for analysis, eliminating the need for sample transfer to individual vials. Collection plate manufacturers provide a variety of options to seal plates. The analyst must prioritize the assay needs (evaporation, contamination, freezing, storage, resealable, etc.) and select the appropriate sealing method.

Reversed Phase Extractions

The small bed mass of sorbent in the disk allows for the use of small solvent volumes. A general guide to solvent volumes for a disk SPE method using reversed phase sorbents (C18, C8, C2 or Universal Resin) is listed in this figure. Each assay will need further optimization in terms of selecting the best wash solvent composition (10% methanol, for instance, will not be optimal for all assays) and the elution solvent. It is recommended that the volume of elution solvent always be optimized, as it will vary depending on the analyte, its affinity for the chosen sorbent, and the strength of the eluting solvent.

**Example SPE Method for Empore™
C18, C8, C2 or Universal Resin
Disk Plates**

<p>Condition 100 µL Methanol 200 µL Water or Buffer</p>
<p>Load 96 Samples into Disk Plate</p>
<p>Wash 500 µL Water 200-500 µL Methanol/Water (10/90)</p>
<p>Elute 100-150 µL Organic Solvent Dilute with Aqueous Component* (commonly 1-3 times volume of organic used)</p>
<p>Direct Inject Aliquot LC/MS/MS</p>

* Resulting mixture is compatible with mobile phase for direct injection

Mixed Phase Extractions

A general solvent volume guide is listed below for a disk SPE method using mixed phase cation (MPC) exchange sorbent. The MPC sorbent is silica bonded with both octyl (C8) and benzene sulfonic acid groups. The conditioning, load and wash volumes should not vary much. However, it is recommended that the volume of elution solvent be optimized to identify the minimum volume that will elute the drug reproducibly from the phase. Two common elution solvent choices are shown. Alternate solvents which may be used include methanol/NH₄OH or acetonitrile/NH₄OH (98/2, v/v).


**Example of SPE Method
Using Empore™ Mixed Phase
Cation Disk Plates**

<p>Condition 100 µL Methanol 200 µL Water 200 µL Buffer</p>
<p>Load 96 Samples into Disk Plate (Typically 50:50 Sample: Buffer, v/v)</p>
<p>Wash 500 µL Water 500 µL Acetic Acid 1.0 M 500 µL Methanol</p>
<p>Elute 150-300 µL ACN: IPA: NH₄OH or CH₂Cl₂:IPA: NH₄ OH (78:20:2, v/v/v)</p>
<p>Evaporate Solvent Reconstitute In Mobile Phase</p>
<p>Inject Aliquot LC/MS/MS</p>

Conditioning Approaches


The thin bed volume of the Empore™ disk allows for the use of small solvent volumes. Given that conditioning can be performed with 4 to 5 times the bed volume, a 100 μL volume of methanol is sufficient to condition the disk. Methanol quickly soaks into the prefilter and disk, and water (200 μL) can be immediately added without a vacuum step between the two additions. Larger volumes of methanol and water can be used, but a vacuum step is necessary after each addition. The disk is reported to be less susceptible to drying out during conditioning than traditional packed particle beds.

An advanced conditioning scheme which further exploits the unique properties of the Empore membrane was published by Pfizer [Janiszewski et al., *Rapid Communications in Mass Spectrometry*, 11, 1033-1037 (1997)]. In the reported assay, a semi-automated liquid handling workstation (Tomtec Quadra 96) allowed the dispensing of methanol (50 μL), followed by an air gap (25 μL), then serum (300 μL) in rapid succession from the same pipette tip. A water step was omitted and vacuum was on during this addition. This novel technique can be evaluated for different analytes, sample matrices, and sorbent types for increased sample prep throughput.



Eliminating the Evaporation Step

The small bed volume of the Empore disk format allows for the use of reduced elution volumes, often eliminating the need for a time-consuming evaporation and reconstitution step. A common approach is to elute with a small volume of organic solvent, then add a volume of aqueous liquid so that the composition of the resulting solution is compatible with mobile phase. Another approach is to elute using a solvent with sufficient organic content to desorb analyte but which is also compatible with mobile phase for direct injection.




Evaporating Solvent Before Reconstitution

Some assays [e.g., mixed phase cation (MPC)] may require the dry-down of eluate (using nitrogen gas, either with or without a heating block) and reconstitution in another solvent suitable for injection. Options are commercially available from a number of suppliers for 96-well plate dry-down.

For more information on enhancing approaches to conditioning, washing, loading and eluting, please refer to the Empore™ High Performance Extraction Disk Plate “Method Optimization Guide.”

Identifying Particle Lots

Each plate has a label with a series of identifying letters and numbers. The first six digits in the identifier contain lot number information. In the example CB003601211-032-1, CB translates to the sorbent type (C18-SD) and 0036 refers to the membrane lot number. The numbers following CB0036, e.g., -01211-032-1, refer to manufacturing process information used to track the raw materials and procedures used to assemble each plate. **Important Note:** Several membrane lots can be made from one particle lot of material, and it is commonly the particle lot that is of concern when identifying enough product to complete an entire analytical study. Contact 3M for further details.



Warranty Information: All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M WARRANTS ONLY THAT ITS PRODUCTS WILL MEET 3M SPECIFICATIONS AT THE TIME OF SHIPMENT. THE FOREGOING WARRANTY IS MADE IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND FREEDOM FROM NON-INFRINGEMENT. **Empore™ Sample Preparation Products are intended for solid phase extraction during scientific research only. These products are not intended or warranted for use in medical devices or in assessment and treatment of clinical patients.**

Limitation of Remedies: If Products are proven not to meet 3M specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: THE REMEDIES PROVIDED HEREIN ARE EXCLUSIVE REMEDIES AGAINST 3M FOR ANY ALLEGED OR ACTUAL NONCONFORMANCE TO SPECIFICATIONS OR DEFECT OR OTHER FAILURE IN PRODUCTS. UNDER NO CIRCUMSTANCES IS 3M LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING LOST PROFITS) IN ANY WAY RELATED TO THE PRODUCT OR SUPPLY OF PRODUCT UNDER ANY THEORY OF LAW INCLUDING, BUT NOT LIMITED TO CONTRACT, NEGLIGENCE AND STRICT LIABILITY.



Bioanalytical Technologies Project

3M Center, Building 270-02A-08
St. Paul, MN 55144-1000
1-888-509-5330
Fax: 1-651-736-4882
Website: www.3M.com/empore