


# BluePippin™ Cassette Kit Reference Chart

Software v.6.20 Cassette Definition Set 18



**For Size Selection:** use **0.60** for optical calibration setting

Product Number Dye-Free**	Size Range <sup>#</sup>	Collection Times	Minimum size Distribution as Expressed by (CV)	Electrophoresis Voltage	Gel Cassette Type	Marker Type	Marker
BHZ7510	50 kb target only	9-10 hr	20%	Pulsed field	0.75% agarose	External	Z1
BMF7510	18 – 27 kb	3.8 – 4.3 hr	20%	Pulsed field	0.75% agarose	External	T1
BUF7510	10 – 18 kb	2.9 – 3.7 hr	20%	Pulsed field	0.75% agarose	External	U1
PAC30KB††	>30 – 40 kb†	6.0 hr	NA	Pulsed field	0.75% agarose	External	U1
PAC20KB††	>4 – 20 kb†	2.0 hr	NA	Pulsed field	0.75% agarose	External	S1
BLF7510	>4 – 20 kb†	2.0 hr	NA	Pulsed field	0.75% agarose	External	S1
BLF7510	3 – 10 kb	2.3 – 3.1 hr	20%	Pulsed field	0.75% agarose	External	S1
BLF7510	2 – 6 kb	1.4 – 2.4 hr	20%	Pulsed field	0.75% agarose	External	S1
BLF7510	1 – 6 kb	3 – 8 hr	10%	25 V, constant	0.75% agarose	External	S1
BDF1510	250 – 1500 bp	28 – 75 min	8%	150 V, constant	1.5% agarose	Internal	R2
BDF2010	100 – 600 bp	50 – 100 min	8%	100 V, constant	2% agarose	Internal	V1
BEF2010	100 – 600 bp	50 – 100 min	8%	100 V, constant	2% agarose	External	M1
BDF3010	90 – 250 bp	50 – 100 min	8%	100 V, constant	3% agarose	Internal	Q2

\*\* Sample DNA is not detected. Fluorescently labeled DNA markers are used to determine collection timing.

# Size range specifications are based on "Tight" mode collections, see next page.

† "High-Pass" protocol. Collects remaining DNA fragments above target size value.

†† "High-Pass" protocols, specifically for PacBio® SMRTbell™ template protocols

# BluePippin™ Cassette Kit Reference Chart

Software v.6.20 Cassette Definition Set 18



Cassette Kit Product Number	Cassette Type Software Definitions	DNA Control Kit Product Number
BMF7510	0.75% DF 18 – 27 kb Marker T1	CON7604
BUF7510	0.75% DF 10 – 18 kb Marker U1	CON7604
BLF7510	0.75% DF 2 – 6 kb Marker S1	CON7504
BLF7510	0.75% DF 3– 10 kb Marker S1	CON7504
BLF7510	0.75% DF Low Voltage 2 – 6 kb Marker S1	CON7504
BDF1510	1.5% DF Marker R2	CIS1504
BDF2010	2% DF Marker V1	CIS2004
BEF2010	2% EF Marker M1	CON2004
BDF3010	3% DF Marker Q2	CIS3004
PAC20KB, BLF7510	High Pass Protocols	CHP7504

	DNA Target Size Ranges	
	Tight Mode	Range Mode
0.75% agarose	40 or 50 kb*	NA
0.75% agarose	18 – 27 kb	13 – 33 kb
0.75% agarose	10 – 18 kb	9 – 20 kb
0.75% agarose	5 – 10 kb	3 – 12 kb
0.75% agarose	2 – 6 kb	1 – 7 kb
1.5% agarose	250 – 1500 bp	200 – 1700 bp
3% agarose	90 – 250 bp	80 – 270 bp

\* Size selection may be made at these target values only. However the selection width may be widened to increase sample yield.



978.922.1832 [www.sagescience.com](http://www.sagescience.com)