

# A REVOLUTION IN HbA1c MEASUREMENT



Affordable A1c  
testing

Semi-automated  
analyser

Easy to use

Results in just  
4 minutes

**EKF**  
**DIAGNOSTICS**

SCIENCE MADE SIMPLE

## EASY TO USE

- User-friendly features minimise training time
- Step by step instructions on display
- User selectable language
- Just 4 µl blood required from finger prick or venous sample
- Innovative blood collector allows easy and consistent sampling

## FAST AND ACCURATE

- Results within 4 minutes
- Uses boronate affinity methodology which is widely recognised as interference-free
- Measuring range: 4-15% A1c DCCT
- Imprecision: CV < 3% at 7% A1c DCCT
- Unaffected by Hb variants, which do not result in reduced erythrocyte life span
- Traceable to the IFCC reference method

## EFFICIENT DATA HANDLING

- Barcode reader to scan calibration data, patient and operator ID
- Stores up to 7,000 results
- User selectable dual reporting. User can select % DCCT, IFCC mmol/mol, eAG mg/dl or eAG mmol/l
- USB port
- Optional label printer available

## ERGONOMIC DESIGN

- Compact: 95 mm(h) x 205 mm(w) x 135 mm(d)
- Lightweight: 0.7 kg



**1**  
Scan lot calibration barcode. Remove foil from cartridge and insert into analyser.



**2**  
When prompted press reagent bead into cartridge using the reverse end of the sampling tool.



**3**  
Insert sample tool into cartridge and snap off handle.



**4**  
Close analyser lid. Result displayed within 4 minutes.

Affordable, accurate and reliable **HbA1c measurement**



### Distributed by

### Manufacturer

**Quotient Diagnostics Ltd**  
Russell House, Molesey Road,  
Walton-on-Thames, Surrey, UK  
KT12 3PE  
**Tel.** +44 (0)1932 220 124



ISO13485:  
2003 Certified.  
Registration ref:  
MD 505902

### Sales

**EKF Diagnostics Holdings plc**  
Avon House, 19 Stanwell Road,  
Penarth, Cardiff, UK CF64 2EZ  
**Tel.** +44 (0)29 20 710 570  
[www.ekfdiagnostics.com](http://www.ekfdiagnostics.com)

**EKF**  
**DIAGNOSTICS**

SCIENCE MADE SIMPLE

Disclaimer: Quotient Diagnostics Ltd, an EKF Diagnostics Holdings plc company, reserves the right to change the product specifications without prior notice.

Revision 3.0-08/2013