

Poloxamer 188, 10% Solution

Product Data Sheet

Product information for MIR 6230

SDS and Certificate of Analysis available at mirusbio.com/6230



SPECIFICATIONS

Storage	Store Poloxamer 188 Solution at 15-30°C.
Product Guarantee	As indicated on the product label, when properly stored and handled.

Product Description:	Poloxamer 188 is a non-ionic detergent that is added to batch mammalian cell cultures to reduce foaming and improve the resistance of cell membranes to hydrodynamic shearing. Poloxamer 188 is 0.1 µm sterile-filtered and is provided at a concentration of 10%. The recommended working concentration of Poloxamer 188 is generally 0.1-0.3% (10-30 mL/L) of culture medium.
Product Size:	100 mL
Form/Appearance:	Clear, colorless solution
Molecular Formula:	(C ₃ H ₆ O.C ₂ H ₄ O) _x
CAS Number:	9003-11-6
Precautions:	Review the product-specific Safety Data Sheet (SDS) for additional information before handling this product.

For Research Use and Further Manufacturing; Not for Administration into Humans

Poloxamer 188, 10% Solution is a filter-sterilized, animal-origin-free reagent that is manufactured under cGMP conditions in an ISO-compliant facility.

► NOTES



SDS and Certificate of Analysis available at [mirusbio.com/6230](https://www.mirusbio.com/6230)

©1996-2020 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners.
For terms and conditions, visit www.mirusbio.com

Rev.A 123119

Mirus Bio LLC

www.mirusbio.com | techsupport@mirusbio.com | U.S. Toll Free: 844.MIRUSBIO (844.647.8724) | Direct: +1.608.441.2852